



**MEDISCA**<sup>®</sup>  
YOUR TRUSTED PARTNER IN COMPOUNDING

**PRODUCT NO.: 6707**



## **RAPID DISSOLVE TABLET MOLD (RDT), PROFILLER**

96 Cavities, 0.940 mL/cavity

**DESCRIPTION:** This device can be used to make 96 rapid dissolve tablets simultaneously and can also be used for commercially available bases.

This mold is dishwasher safe when cleaned with water only.

### **TECHNICAL SPECIFICATIONS:**

Overall Dimensions, W x D x H (mm / in)	235 x 330 x 205 / 9.3 x 13.0 x 0.8	
Weight, Tamper and Cavity Plate (kg / lbs)	3.35 / 7.38	
Individual Cavity, Ø x H (mm / in)	13.4 x 7.5 / 0.53 x 0.30	
Hole volume (mL)	0.940	
Material	Tamper Plate	Hard anodized aluminum
	Cavity Plate	Teflon coated Hard anodized aluminum